Serum CrossLaps® ELISA

For the quantification of degradation products of C-terminal telopeptides of Type I collagen in human serum and plasma

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IDS Ltd is not responsible for any other use of the kit or consequence hereof than the one specified above. Neither for misuse, e.g. use deviating from the procedure described in this manual. Furthermore, IDS Ltd is not to be made responsible for any diagnoses or conclusions made by the user or third party based on the results obtained with the Serum CrossLaps® ELISA kit nor for any consequences such interpretations may cause.

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INTRODUCTION

Intended use

The Serum CrossLaps® ELISA is an enzyme immunological test for the quantification of degradation products of C-terminal telopeptides of Type I collagen in human serum and plasma.

The Serum CrossLaps® ELISA assay is intended for *in vitro* diagnostic use as an indication of human bone resorption and may be used as an aid in

A. Monitoring bone resorption changes of

- 1) Anti resorptive therapies in postmenopausal women:
 - a) Hormone Replacement Therapies (HRT) with hormones and hormone like drugs
 - b) Bisphosphonate therapies
- 2) Anti resorptive therapies in individuals diagnosed with osteopenia;
 - a) Hormone Replacement Therapies (HRT) with hormones and hormone like drugs
 - b) Bisphosphonate therapies

B. Predicting skeletal Response (Bone Mineral Density) in postmenopausal women undergoing anti resorptive therapies

- a) Hormone Replacement Therapies (HRT) with hormones and hormone like drugs
- b) Bisphosphonate therapies

Limitations

The use of the test has not been established to predict the development of osteoporosis or future fracture risk.

The use of the test has not been established in hyperparathyroidism or hyperthyroidism.

When using the test to monitor therapy, results may be confounded in patients afflicted with clinical conditions known to affect bone resorption e.g. bone metastases, hyperparathyroidism or hyperthyroidism.

Serum CrossLaps® ELISA results should be interpreted in conjunction with clinical findings and other diagnostic results and should not be used as a sole determinant in initiating or changing therapy. Do not interchange Serum CrossLaps® ELISA values with Urine CrossLaps® ELISA values.

Summary and explanation of the test

Type I collagen accounts for more than 90% of the organic matrix of bone and is synthesized primarily in bone (1). During renewal of the skeleton, Type I collagen is degraded, and small peptide fragments are excreted into the bloodstream. These fragments can be measured by Serum CrossLaps® ELISA. The measurements of the specific degradation products of Type I collagen in both urine (2) and serum (3) by a competitive CrossLaps have been reported.

The sandwich assay has been reported as useful for follow up of anti resorptive treatment of patients with metabolic bone diseases (3-17).

Principle of the procedure

The Serum CrossLaps® ELISA is based on two highly specific monoclonal antibodies against the amino acid sequence of EKAHD- β -GGR, where the aspartic acid residue (D) is β -isomerized. In order to obtain a specific signal in the Serum CrossLaps® ELISA, two chains of EKAHD- β -GGR must be cross linked. Standards, control, or unknown serum samples are pipetted into the appropriate microtitre wells coated with streptavidin, followed by application of a mixture of a biotinylated antibody and a peroxidase conjugated antibody. Then, a complex between CrossLaps antigens, biotinylated antibody and peroxidase conjugated antibody is generated, and this complex binds to the streptavidin surface via the biotinylated antibody. Following the one step incubation at room temperature, the wells are emptied and washed. A chromogenic substrate is added and the colour reaction is stopped with sulfuric acid. Finally, the absorbance is measured.

PRECAUTIONS

The following precautions should be observed in the laboratory:

- Do not eat, drink, smoke or apply cosmetics where immunodiagnostic materials are being handled
- Do not pipette by mouth.
- Wear gloves when handling immunodiagnostic materials and wash hands thoroughly afterwards
- Cover working area with disposable absorbent paper

Warnings

For *in vitro* use only.

- All reagents and laboratory equipment should be handled and disposed of as if they were infectious.
- Do not use kit components beyond the expiry date and do not mix reagents from different lots.

HAMA interference

Some individuals have antibodies to mouse immunoglobulins (HAMA), which can cause interference in immunoassays that employ murine monoclonal antibodies, such as Serum CrossLaps. In rare cases, the content of HAMA exceeds the capacity of the blocking agent incorporated into Serum CrossLaps leading to a false-positive test result. Therefore, Serum CrossLaps values should be used only in conjunction with information available from the clinical evaluation of the patient.

Storage

Store the Serum CrossLaps® ELISA kit upon receipt at 2-8°C. Under these conditions the kit is stable up to the expiry date stated on the box.

MATERIALS

Specimen collection

Collect blood by venipuncture taking care to avoid haemolysis. Separate the serum from the cells within 3 hours after collection of blood. It is recommended to freeze (<-18°C) samples immediately. For optimal results it is recommended to draw blood as fasting morning samples (18). Also for monitoring the individual patient, follow up samples should be collected under same conditions as the baseline sample.

When analysing plasma, both heparia and EDTA plasma may be used.

Materials supplied

Before opening the kit, read the section on Precautions. The kit contains reagents sufficient for 96 determinations.

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Microwell strips (12x8 wells) pre-coated with streptavidin. Supplied in a plastic frame.

CrossLaps Standard **CAL 0**

One vial (min. 5.0 mL/vial) of ready for use PBS buffered solution with protein stabiliser and preservative.

CrossLaps Standards **CAL 1-5**

Five vials (min. 0.4 mL/vial) of ready for use, CrossLaps standard in a PBS buffered solution with protein stabiliser and preservative. The exact value of each Standard is printed on the QC Report.

Control **CTRL 1-2**

Two vials (min. 0.4 mL/vial) of ready for use, desalted urinary antigens of human origin in a PBS buffered solution with protein stabiliser and preservative. Please refer to enclosed QC Report for control range.

Biotinylated Antibody **Ab BIOTIN**

One vial (min. 0.25 mL) of a concentrated solution of a biotinylated monoclonal murine antibody specific for degradation products of C-terminal telopeptides of Type I collagen, raised. Prepared in a buffered solution with protein stabiliser and preservative.

Peroxidase Conjugated Antibody [ENZYMCONJ]

One vial (min. 0.25 mL) of a concentrated solution of a peroxidase conjugated murine monoclonal antibody specific for degradation products of C-terminal telopeptides of Type I collagen. Prepared in a buffered solution with protein stabiliser and preservative.

Incubation Buffer [BUF]

One vial (min. 19 mL) of a ready for use buffered solution with protein stabiliser, detergent and preservative.

Substrate Solution [SUBS] [TMB]

One vial (min. 12 mL) of a ready for use tetramethylbenzidine (TMB) substrate in an acidic buffer. Please note that the chromogenic substrate might appear slightly blueish.

Stopping Solution [H2SO4]

One vial (min. 12 mL) of ready for use 0.18 mol/L sulfuric acid.

Washing Buffer [WASHBUF] [50x]

One vial (min. 20 mL) of a concentrated washing buffer with detergent and preservative.

Sealing tape

Adhesive film for covering wells during incubation.

Materials required — not supplied

- Containers for preparing the Antibody Solution and the Washing Solution
- Precision micropipettes to deliver 50-200 µL
- Distilled water
- Precision 8 or 12 channel multipipette to deliver 100 µL, and 150 µL
- Microwell mixing apparatus
- Microtiter plate reader

ASSAY PROCEDURE

Mix all reagents and samples before use (avoid foam)

Determine the number of strips needed for the assay. It is recommended to test all samples in duplicate. In addition, for each run a total of 16 wells are needed for the standards and controls. Place the appropriate number of strips in the plastic frame. Store unused immuno strips in the tightly closed foil bag with desiccant capsules.

Prior to use, prepare and equilibrate all solutions to room temperature. **Perform the assay at room temperature (18-22°C).**

1 Preparation of the Antibody Solution:

ATTENTION: Prepare the following Antibody Solution maximum 30 minutes before starting the assay. Mix the Biotinylated Antibody [Ab BIOTIN], Peroxidase Conjugated Antibody [ENZYMCONJ] and Incubation Buffer [BUF] in the volumetric ratio 1+1+100 in an empty container. Mix carefully and avoid formation of foam. **Prepare a fresh solution before each run of the assay.**

2 One Step incubation

Pipette 50 µL of either **Standards** [CAL 0]-[5], **Control** [CTRL 1]-[2], or unknown samples into appropriate wells followed by 150 µL of the **Antibody Solution**. Cover the immunostrips with sealing tape and incubate for 120±5 minutes at room temperature (18-22°C) on a microtitre plate mixing apparatus (300 rpm).

3 Washing

Wash the immunostrips 5 times manually with 300 µL diluted **Washing Buffer** ([WASHBUF] [50x] diluted 1+50 in distilled water). Using an automated plate washer, follow the instructions of the manufacturer or the guidelines of the laboratory. Usually 5 washing cycles are adequate. Make sure that the wells are **completely emptied** after each manual or automatic washing cycle.

4 Incubation with chromogenic substrate solution

Pipette 100 µL of the **Substrate Solution** [SUBS] [TMB] into each well and incubate for 15±2 minutes at room temperature (18-22°C) in the dark on the mixing apparatus (300 rpm). Use sealing tape. Do not pipette directly from the vial containing TMB substrate but transfer the needed volume to a clean reservoir. Remaining substrate in the reservoir should be discarded and not returned to vial TMB.

5 Stopping of colour reaction

Pipette 100 µL of the **Stopping Solution** [H2S04] into each well.

6 Measurement of absorbance

Measure the absorbance at 450 nm with 650 nm as reference within two hours.

Limitations of the procedure

If the absorbance of a sample exceeds that of **Standard 5**, the sample should be diluted in **Standard 0** and re-analysed.

QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

RESULTS

Calculation of results

A quadratic curve fit can be used.

Alternatively, calculate the mean of the duplicate absorbance determinations. Construct a standard curve on graph paper by plotting the mean absorbances of the six standards 0-5 (ordinate) against the corresponding CrossLaps concentrations (abscissa). Determine the CrossLaps concentration of the controls and each patient sample by interpolation.

Example of results obtained:

Standards/ Controls/ Samples	CrossLaps conc. (ng/mL)	A ₄₅₀₋₆₅₀ (nm) Obs 1/ Obs 2	Mean A ₄₅₀₋₆₅₀ (nm)	Interpolated CrossLaps conc. (ng/mL)
Standard 0	0.000	0.066 / 0.065	0.066	
Standard 1	0.178	0.210 / 0.209	0.210	
Standard 2	0.489	0.472 / 0.448	0.460	
Standard 3	0.960	0.844 / 0.819	0.832	
Standard 4	1.902	1.598 / 1.560	1.579	
Standard 5	2.494	2.061 / 2.004	2.033	
Control 1		0.349 / 0.354	0.352	0.355
Control 2		0.918 / 0.952	0.935	1.086
Sample I		0.140 / 0.138	0.139	0.091
Sample II		0.447 / 0.439	0.443	0.469
Sample III		1.305 / 1.303	1.304	1.555

Please note: The data above are for illustration only and should not be used to calculate the results of any run.

Performance characteristics

Detection limit: 0.020 ng/mL CrossLaps

This is the concentration corresponding to three standard deviations above the mean of 21 determinations of the blank (CrossLaps Standard 0).

Imprecision

The imprecision of the Serum CrossLaps® ELISA was evaluated for three serum samples. The results are summarised in the table below:

InterAssay Variation (n=10)

Mean (ng/mL)	SD (ng/mL)	CV (%)
0.121	0.013	10.9
0.444	0.043	9.7
1.967	0.050	2.5

IntraAssay Variation (n=10)

Mean (ng/mL)	SD (ng/mL)	CV (%)
0.121	0.004	3.0
0.444	0.007	1.7
1.967	0.035	1.8

Dilution/Linearity

The Serum CrossLaps® ELISA is linear in the range 0.020 ng/mL to 3.380 ng/mL of CrossLaps.

Serum samples with the concentration of 0.460 - 0.668 ng/mL CrossLaps were diluted with standard 0 and the concentration of CrossLaps were determined with Serum CrossLaps® ELISA. The serum neat sample is set to 100%.

The data below is calculated from 3 different runs:

Dilution Procedure		
Serum [%]	Standard 0 [%]	Recovery [% of expected value]
100	0	100
90	10	103
80	20	103
70	30	101
60	40	102
50	50	105
40	60	109
30	70	107
20	80	97
10	90	100
Mean		103

Interference:

The interference of Ditaurobilirubin, Hemoglobin and IntraLipid on the measurement of CrossLaps in serum by Serum CrossLaps® ELISA was investigated.

In the concentration listed below no interference was detected:

- Ditaurobilirubin: 600 mg/L
- Hemoglobin: 10 g/L
- IntraLipid: 10 g/L

Expected values

It is advisable for a laboratory to establish its own range of normal and pathological values. As an example, the mean values and standard deviations for various populations are given below. For further reading, please refer to the reference list. All samples were morning fasting samples from healthy individuals.

Populations	Number of subjects	Mean Values* (ng/mL)	95% Confidence Interval (ng/mL)
Post-menopausal women	193	0.439	0.142 – 1.351
Pre-menopausal women	226	0.287	0.112 – 0.738
Males	125	0.294	0.115 – 0.748

Day to Day Individual Variation

The Day to Day Intra individual Variation was assessed by analyzing serum samples (morning fasting) from 11 healthy post menopausal women at five time points over 2 weeks

Subject No	Visit 1	Visit 2	Visit 3	Visit 4	Visit 5	Mean (ng/mL)	SD (ng/mL)	CV (%)
1	0,423	0,428	0,396	0,460	0,445	0,430	0,024	6
2	0,461	0,523	0,500	0,535	0,539	0,512	0,032	6
3	0,850	0,731	0,761	0,782	0,764	0,778	0,045	6
4	0,377	0,468	0,455	0,499	0,440	0,448	0,045	10
5	0,918	0,834	0,791	0,781	0,714	0,808	0,075	9
6	0,268	0,249	0,246	0,257	0,258	0,255	0,009	3
7	0,431	0,457	0,468	0,494	0,506	0,471	0,030	6
8	0,666	0,587	0,670	0,595	0,728	0,655	0,063	10
9	0,323	0,357	0,341	0,409	0,345	0,355	0,033	9
10	0,419	0,520	0,541	0,491	0,470	0,488	0,047	10
11	0,353	0,472	0,429	0,464	0,400	0,424	0,049	12

CLINICAL DATA

The Serum CrossLaps® ELISA has been used to monitor treatment in several clinical studies and the CrossLaps values have been compared to Bone Mineral Density (BMD-spine) measurements.

All the clinical studies presented below were performed according to the European Standard for good clinical practice (GCP and GLP).

Most of the clinical studies presented here were conducted on white Danish women. However, several studies have been published showing that other demographic groups display similar CrossLaps decrease in response to anti resorptive therapies (9-13).

For all the data presented below fasting morning samples have been used.

The Bone mineral density was measured at the Lumbar spine (L1, L4).

The change in the bone mineral density is presented below α -BMD. α -BMD is defined as the slope of the linear regression line for BMD-spine versus time (years) for the period of treatment. In most cases the calculation of α BMD involves a minimum of 5 BMD-spine measurements. The α -BMD thus represents the % change in BMD-spine per year.

Because there to this date is no universal agreement as to what constitutes positive BMD response we have calculated the sensitivity and specificity using two different cut off values for α -BMD; α -BMD>0 and α -BMD>1.

The sensitivity is defined as the percent of the study population with a positive BMD response and who have a % change from baseline of Serum CrossLaps® ELISA which is 40% or greater.

The specificity is defined as the percent of the study population without a positive BMD response and who have a % change from baseline of Serum CrossLaps® ELISA that is less than 40%.

Bisphosphonate studies

Below is shown the Serum CrossLaps® ELISA data from two different bisphosphonate studies.

Alendronate

- Women between age 40 and 59 years, 6 months to 3 years since menopause
- 12 participants on placebo (500 mg calcium)
- 42 participants on active treatment (5 mg (n=16), 10 mg (n=14), 20 mg (n=12)) Alendronate and 500 mg calcium)
- Treatment period: 2 - 3 years

Serum CrossLaps® ELISA


	Placebo group			Alendronate group		
	Mean (ng/mL) (95% Confidence Int.)	SD	SEM	Mean (ng/mL) (95% Confidence Int.)	SD	SEM
Baseline	0.676 (0.543-0.809)	0.235	0.068	0.603 (0.540-0.666)	0.208	0.032
After 6 months treatment	0.623 (0.498-0.748)	0.221	0.064	0.191 (0.144-0.238)	0.153	0.024

Ibandronate

- Women less than 75 years, more than ten years after menopause and have a BMD forearm 1.5 SD or more below the standard for healthy pre menopausal women
- 17 participants on placebo (1000 mg calcium)
- 36 participants on active treatment: (2.5 mg (n=20), 5 mg (n=16)) ibandronate and 1000 mg calcium
- Treatment period: 1 year

Serum CrossLaps® ELISA

	Placebo group			Ibandronate group		
	Mean (ng/mL) (95% Confidence Int.)	SD	SEM	Mean (ng/mL) (95% Confidence Int.)	SD	SEM
Baseline	0.590 (0.502-0.678)	0.185	0.045	0.614 (0.536-0.692)	0.240	0.040
After 6 months treatment	0.325 (0.241-0.409)	0.178	0.043	0.136 (0.093-0.179)	0.131	0.022



Using a cut-off for Serum CrossLaps® ELISA of 40% change from baseline the following sensitivities, specificities and 95% confidence intervals are obtained.

Below is indicated the actual fractions:

	Ibandronate		Alendronate	
	Sensitivity	Specificity	Sensitivity	Specificity
α -BMD>0	93% (81-99) 40/43	30% (7-65) 3/10	89% (73-97) 31/35	92% (64-100) 12/13
α -BMD>1	94% (81-99) 33/35	22% (6-48) 4/18	90% (73-98) 26/29	68% (43-87) 13/19

Below is shown distribution plots for the combined bisphosphonate studies. The number over each bar indicates the number of participants in each class.

HRT Studies

Below is shown the Serum CrossLaps® ELISA data from three different HRT studies.

Tibolone

- Women less than 75 years and more than ten years after meno pause
- 13 participants on placebo (400 mg calcium/day)
- 49 participants on active treatment (1.25 mg (n=25) or 2.5 mg (n=24) Tibolone and 400 mg calcium/day
- Treatment period 2 years

Serum CrossLaps® ELISA

	Placebo group			Tibolone group		
	Mean (ng/mL) (95% Confidence Int.)	SD	SEM	Mean (ng/mL) (95% Confidence Int.)	SD	SEM
Baseline	0.264 (0.217-0.311)	0.085	0.024	0.339 (0.302-0.376)	0.130	0.019
After 6 months treatment	0.287 (0.232-0.342)	0.099	0.028	0.192 (0.161-0.223)	0.113	0.016

HRT I

- Women more than 45 years, 1 to 6 years since menopause
- 42 participants on placebo (400 mg calcium/day)
- 120 participants on active treatment:

Days 1-16	Days 17-28
E 1 mg	E 1 mg + G 25 µg
E 2 mg	E 2 mg + G 25 µg
E 2 mg	E 2 mg + G 50 µg
E 1 mg + G 25 µg	E 1 mg + G 25 µg continuously
- E = estradiol 17β, G = gestodene, active treatment also receive 400 mg calcium/day
- Treatment period: 2 years.

Serum CrossLaps® ELISA


	Placebo group			HRT I group		
	Mean (ng/mL) (95% Confidence Int.)	SD	SEM	Mean (ng/mL) (95% Confidence Int.)	SD	SEM
Baseline	0.389(0.348-0.430)	0.136	0.021	0.411(0.386-0.437)	0.140	0.013
After 6 months treatment	0.396(0.349-0.443)	0.153	0.024	0.182(0.164-0.200)	0.098	0.009

HRT II

- Women between 65 and 70 years and BMC forearm below 1 SD of healthy pre-menopausal women
- 17 participants on placebo (1000 mg calcium/day)
- 15 participants on active treatment: 50 µg estradiol, 1 mg norethisterone and 1000 mg calcium/day

Serum CrossLaps® ELISA


	Placebo group			HRT II group		
	Mean (ng/mL) (95% Confidence Int.)	SD	SEM	Mean (ng/mL) (95% Confidence Int.)	SD	SEM
Baseline	0.350(0.301-0.399)	0.099	0.025	0.371(0.306-0.436)	0.135	0.033
After 6 months treatment	0.293(0.240-0.346)	0.106	0.027	0.159(0.108-0.210)	0.106	0.026



Using a cut off for Serum CrossLaps® ELISA of 40% change from baseline the following sensitivities, specificities and 95% confidence intervals are obtained. Below is indicated the actual fractions:

	Tibolone		HRT I		HRT II	
	Sensitivity	Specificity	Sensitivity	Specificity	Sensitivity	Specificity
α -BMD>0	58% (44-72) 31/53	90% (56-100) 9/10	80% (72-87) 16/20	90% (77-97) 37/41	60% (39-79) 15/25	100% (59-100) 7/7
α -BMD>1	65% (49-79) 28/43	80% (56-94) 16/20	83% (74-94) 90/109	78% (64-88) 42/54	62% (38-82) 13/21	91% (59-100) 10/11

Below is shown distribution plots for the combined HRT studies. The number over each bar indicates the number of participants in each class.

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REFERENCES

1. Bjarnason NH. & Christiansen, C. Early Response in Biochemical Markers Predicts Long-term Response in Bone Mass During Hormone Replacement Therapy in Early Postmenopausal Women. *Elsevier* 26; 561-569 (2000).
2. Bonde M, Garnero P, Fledelius C, Qvist P, Delmas PD, Christiansen C. Measurement of bone degradation products in serum using antibodies reactive with an isomerized form of an 8 amino acid sequence of the C-telopeptide of Type I collagen. *J Bone Miner Res* 12:1028-1034 (1997).
3. Bonde M, Qvist P, Fledelius C, Riis BJ, Christiansen C. Applications of an enzyme immunoassay for a new marker of bone resorption (CrossLaps®) - follow up on hormone replacement therapy and osteoporosis risk assessment. *J Clin Endocrinol & Metab* 80:864-8 (1995).
4. Buclin T, Rochat MC, Burckhardt P, Azria M, Attinger M. Bioavailability and biological efficacy of a new oral formulation of salmon calcitonin in healthy volunteers. *J Bone Miner Res* 17;8-1478 (2002).
5. Burgeson RE. New collagens, new concepts. *Annu Rev Cell Biol* 4:551-77 (1988).
6. Chailurkit LO, Aunphongpuwanart S, Ongphiphadhanakul B, Jongjaroenprasert W, Sae-Tung S, and Rajatanavin R. Efficacy of intermittent low dose alendronate in Thai postmenopausal osteoporosis. *Endocr Res* 30:29-36 (2004).
7. Christgau S, Bitsch-Jensen O, Bjarnason NH, Gamwell Henriksen E, Qvist P, Alexandersen P, Bang Henriksen D. Serum CrossLaps for Monitoring the Response in Individuals Undergoing Antiresorptive Therapy. *Elsevier* 26; 505-511 (2000).
8. Christgau S, Rosenquist C, Alexandersen P, Bjarnason NH, Ravn P, Fledelius C, Herling C, Qvist P, Christiansen C. Clinical evaluation of the Serum CrossLaps® ELISA, a new assay measuring the serum concentration of bone-derived degradation products of type I collagen Ctelopeptides. *Clin Chem* 44:11;2290-2300 (1998).
9. Christiansen C, Tanko LB, Warming L, Moelgaard A, Christgau S, Qvist P, Baumann M, Wiczorek L, and Hoyle N. Dose dependent effects on bone resorption and formation of intermittently administered intravenous ibandronate. *Osteoporos Int* (2003).
10. Delmas PD. Markers of bone turnover for monitoring treatment of osteoporosis with antiresorptive drugs. *Osteoporos Int* 11;S66-S76 (2000).
11. Delmas PD, Hardy P, Dain MP. Monitoring Individual Response to Hormone Replacement Therapy with Bone Markers. *Elsevier* 26;553-560 (2000).
12. Fink E, Cormier C, Steinmetz P, Kindermans C, Bouc Le Y, Souberbielle JC. Differences in the Capacity of Several Biochemical Bone Markers to Assess High Bone Turnover in Early Menopause and Response to Alendronate Therapy. *Osteoporos Int* 11;295-303 (2000).
13. Garnero P, Gineyts E, Riou JP, Delmas PD. Assessment of bone resorption with a new marker of collagen degradation in patients with metabolic bone disease. *J Clin Endocrinol & Metab* 79:780-5 (1994).
14. Pedersen BJ, Ravn P, Bonde M. Type I Collagen C-telopeptide Degradation products as Bone Resorption Markers. *J Clin Ligand Assay* 21;2:118-127 (1998).
15. Qvist P, Munk M, Hoyle N, Christiansen C. Serum and plasma fragments of C-telopeptides of type I collagen (CTX) are stable during storage at low temperatures for 3 years. *Clin Chim Acta* 350:167-173 (2004).
16. Ravn P, Neugebauer G, Christiansen C. Association between pharmacokinetics of oral ibandronate and clinical response in bone mass and bone turnover in women with postmenopausal osteoporosis. *Bone* 30;320-324 (2002).
17. Ravn P. & Clemmesen B. Biochemical Markers Can Predict the Response in Bone Mass During Alendronate Treatment in Early Postmenopausal Women. *Bone* 24, 237-244 (1997).
18. Reginster JY, Wilson KM, Dumont E, Bonvoisin B, Barrett J. Monthly Oral Ibandronate is Well Tolerated and Efficacious in Postmenopausal Women: Results From the Monthly Oral Pilot Study. *J Clin Endocrinol Metab* (2005).
19. Reginster JY, Henrotin Y, Christiansen C, Gamwell-Henriksen E, Bruyere Collette J, Christgau S. Bone resorption in post-menopausal women with normal and low BMD assessed with biochemical markers specific for telopeptide derived degradation products of collagen type I. *Calcif Tissue Int* 69;130-137 (2001).
20. Rosen HN, Parker RA, Greenspan SL, Iloputaife ID, Bookman L, Chapin D, Perlmutter I, Kessel B, Qvist P, Rosenblatt M. Evaluation of Ability of Biochemical Markers of Bone Turnover to Predict a Response to Increased Doses of HRT. *Calcif Tissue Int* (2004).
21. Rosen HN, Moses AC, Garber J, Iloputaife ID, Ross DS, Lee SL, Greenspan SL. Serum CTX: A New Marker of Bone Resorption That Shows Treatment Effect More Often Than Other Markers Because of Low coefficient of Variability and Large Changes with Bisphosphonate Therapy. *Calcif Tissue Int* 66;100-103 (2000).

22. Rosenquist C, Fledelius C, Christgau S, Pedersen BJ, Bonde M, Qvist P, Christiansen C. Serum CrossLaps® ELISA. First application of monoclonal antibodies for measurement in serum of bonerelated degradation products from C-terminal telopeptides of type I collagen. *Clin Chem* 44;11:2281-2289 (1998).
23. Tanko LB, Bagger YZ, Alexandersen P, Devogelaer JP, Reginster JY, Chick R, Olson M, Benmamar H, Mindeholm L, Azria M, Christiansen C. Safety and Efficacy of a Novel Salmon Calcitonin (sCT) Technology-Based Oral Formulation in Healthy Postmenopausal Women: Acute and 3-Month Effects on Biomarkers of Bone Turnover. *J Bone Miner Res* 19:1531-1538 (2004).
24. Warming L, Ravn P, Spielman D, Delmas P, Christiansen C. Trimegestone in a low-dose, continuous-combined hormone therapy regimen prevents bone loss in osteopenic postmenopausal women. *Menopause* 11:337-342 (2004).
25. Warming L, Ravn P, Nielsen T, Christiansen C. Safety and efficacy of drospirenone used in a continuous combination with 17beta-estradiol for prevention of postmenopausal osteoporosis. *Cimacteric* 7:203-111 (2004).
26. Warming L, Ravn P, Christiansen C. Levonorgestrel and 17beta-estradiol given transdermally for the prevention of postmenopausal osteoporosis. *Maturitas* 50:78-85 (2005).

Doc: AC-02PL-A
Issue: 5
Date:27 April 2010

 EXP	GB <i>Use By</i> DE <i>Verwendbar bis</i> ES <i>Fecha de caducidad</i> IT <i>Utilizzare entro</i> FR <i>Utiliser jusque</i> NL <i>Houdbaar tot</i> DK <i>Holdbar til</i> CZ <i>Použitelné do</i> SK <i>Použitelné do</i> GR <i>Ημερομηνία λήξης</i> PT <i>Prazo de validade</i> HU <i>Felhasználható</i> SE <i>Använd före</i> PL <i>Użyć przed</i>	 LOT	GB <i>Batch code</i> DE <i>Chargenbezeichnung</i> ES <i>Código de lote</i> IT <i>Codice del lotto</i> FR <i>Code du lot</i> NL <i>Lot nummer</i> DK <i>Lotnummer</i> CZ <i>Číslo šarže</i> SK <i>Číslo šarže</i> GR <i>Αριθμός Παρτίδας</i> PT <i>Código do lote</i> HU <i>Sarzszzám</i> SE <i>Lot nummer</i> PL <i>Kod partii</i>
 REF	GB <i>Catalogue number</i> DE <i>Bestellnummer</i> ES <i>Número de catálogo</i> IT <i>Numero di catalogo</i> FR <i>Référence du catalogue</i> NL <i>Catalogus nummer</i> DK <i>Katalognummer</i> CZ <i>Katalogové číslo</i> SK <i>Katalógové číslo</i> GR <i>Αριθμός καταλόγου</i> PT <i>Referência de catálogo</i> HU <i>Katalógusszám</i> SE <i>Katalognummer</i> PL <i>Numer katalogowy</i>		GB <i>Manufacturer</i> DE <i>Hersteller</i> ES <i>Fabricante</i> IT <i>Fabbricante</i> FR <i>Fabricant</i> NL <i>Fabrikant</i> DK <i>Producent</i> CZ <i>Výrobce</i> SK <i>Výrobca</i> GR <i>Κατασκευαστής</i> PT <i>Fabricante</i> HU <i>Gyártó</i> SE <i>Tillverkare</i> PL <i>Producent</i>
	GB <i>Contains sufficient for <n> tests</i> DE <i>Inhalt ausreichend für <n> Prüfungen</i> ES <i>Contenido suficiente para <n> ensayos</i> IT <i>Contenuto sufficiente per "n" saggi</i> FR <i>Contenu suffisant pour "n" tests</i> NL <i>Inhoud voldoende voor "n" testen</i> DK <i>Indeholder tilstrækkeligt til "n" test</i> CZ <i>Lze použít pro <n> testů</i> SK <i>Obsah postačuje na <n> stanovení</i> GR <i>Περιεχόμενο επαρκές για «n» εξετάσεις</i> PT <i>Conteúdo suficiente para "n" ensaios</i> HU <i>A doboz tartalma <n> vizsgálat elvégzéséhez elegendő</i> SE <i>Räcker till "n" antal tester</i> PL <i>Wystarczy na wykonanie <n> testów</i>	 IVD	GB <i>In Vitro Diagnostic Medical Device</i> DE <i>In-Vitro-Diagnostikum</i> ES <i>Producto sanitario para diagnóstico in vitro</i> IT <i>Dispositivo medico-diagnostico in vitro</i> FR <i>Dispositif médical de diagnostic in vitro</i> NL <i>Medisch hulpmiddel voor in-vitro diagnostiek</i> DK <i>Medicinsk udstyr til in vitro-diagnostik</i> CZ <i>In Vitro diagnostický zdravotnický prostředek</i> SK <i>Zdravotnícka pomocka in vitro</i> GR <i>In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν</i> PT <i>Dispositivo médico para diagnóstico in vitro</i> HU <i>In vitro diagnosztikum</i> SE <i>Medicintekniska produkter för in vitro diagnostik</i> PL <i>Wyrób do diagnostyki In Vitro</i>
	GB <i>Temperature limitation</i> DE <i>Temperaturbegrenzung</i> ES <i>Límite de temperatura</i> IT <i>Limiti di temperatura</i> FR <i>Limites de température</i> NL <i>Temperatuurlimiet</i> DK <i>Temperaturbegrænsning</i> CZ <i>Teplotní rozmezí od do</i> SK <i>Teplotné rozmedzie od do</i> GR <i>Περιορισμοί θερμοκρασίας</i> PT <i>Limites de temperatura</i> HU <i>Hőmérséklettartomány</i> SE <i>Temperaturbegränsning</i> PL <i>Przestrzeżać zakresu temperatury</i>		GB <i>Consult Instructions for Use</i> DE <i>Gebrauchsanweisung beachten</i> ES <i>Consulte las instrucciones de uso</i> IT <i>Consultare le istruzioni per l'uso</i> FR <i>Consulter les instructions d'utilisation</i> NL <i>Raadpleeg de gebruiksaanwijzing</i> DK <i>Se brugsanvisning</i> CZ <i>Viz návod k použití</i> SK <i>Vid' návod na použitie</i> GR <i>Συμβουλευτείτε τις οδηγίες χρήσης</i> PT <i>Consulte as instruções de utilização</i> HU <i>Nézze meg a Használati utasítást</i> SE <i>Se handhavandebeskrivningen</i> PL <i>Sprawdź w instrukcji obsługi</i>

Immunodiagnostic Systems Ltd (IDS Ltd).

UK Immunodiagnostic Systems Ltd (IDS Ltd), 10 Didcot Way, Boldon Business Park, Boldon, Tyne & Wear, NE35 9PD

Tel: +44 (0) 191 519 0660 • Fax: +44 (0) 191 519 0760 • e-mail: info.uk@idsplc.com • www.idsplc.com

USA Immunodiagnostic Systems Inc (IDS Inc.), P.O. Box 17063, Fountain Hills, AZ 85269-7063

Tel: 480-836-7435 • Fax: 480-836-7437 • e-mail: info.us@idsplc.com • www.idsplc.com

Germany Immunodiagnostic Systems GmbH (IDS GmbH), Mainzer Landstrasse 49, 60329 Frankfurt am Main

Tel: +49 (0) 69 3085-5025 • Fax: +49 (0) 69 3085-5125 • e-mail: info.de@idsplc.com • www.idsplc.com

France Immunodiagnostic Systems France SA (IDS France SA) 153 Avenue D'Italie, 75013 PARIS

Tel: +33 (0)1 40 77 04 50 • Fax: +33 (0)1 40 77 04 55 • e-mail: info.fr@idsplc.com • www.idsplc.com

Scandinavia Immunodiagnostic Systems Nordic a/s (IDS Nordic a/s), Marielundvej 30, 2. Sal, 2730 Herlev, Denmark

Tel: +45 44 84 0091 • Fax: +45 44 84 0092 • email: info.nordic@idsplc.com • www.idsplc.com